



ENLIGHTEN™ 4.1.0 Release Notes

August 22, 2024

Introduction

The last formal customer release of ENLIGHTEN™ for Windows, version 3.2.36, was released to the Wasatch Photonics website Oct 12, 2022. The product then went into an extended internal beta for a major refactor and incorporation of a number of new features and capabilities. As a result, the new version 4.1.0 ended up taking longer to complete than anticipated, and also contains a larger-than-usual set of changes and improvements.

This document is provided as a consolidated listing of key changes between 3.2.36 and 4.1.0.

Future releases are expected to be rolled-out at a more regular interval of months rather than years.

Key Changes

The following are high-level changes categorized by change type. For a chronological list of detailed changes, see Appendices at document end.

GUI / Navigation

- GUI now organized into three (3) "Modes" (Raman, Non-Raman and Expert)
 - **Raman** Mode only exposes controls recommended for typical Raman measurements
 - **Non-Raman** Mode exposes controls appropriate for typical Absorbance, Reflectance and other non-Raman measurements
 - **Expert** Mode exposes all controls, including advanced features
- All three modes are organized into five (5) "Views" (Scope, Settings, Hardware, Log, and Factory)
 - **Scope** mode is centered around the large central graph and is intended for real-time data collection
 - **Settings** allows configuration of various application-wide features
 - **Hardware** provides detailed information about connected spectrometer(s), including EEPROM and firmware data
 - **Log** provides a view onto the "enlighten.log" file written to the user's EnlightenSpectra folder
 - **Factory** provides advanced manufacturing data for unit alignment and calibration, and requires authentication via a password
- Notification Bar
 - Always visible
 - Errors highlighted in red
- Presets
 - Users can save common sets of acquisition parameters (integration time, gain, laser power, scan averaging etc) to named "presets" for quick selection when changing samples

- Keyboard Shortcuts
 - Ctrl-E to Edit last-saved thumbnail label
 - Ctrl-T to change integration Time
 - Ctrl-G to change Gain (on supported devices in Expert mode)
 - Ctrl-N to change the Note (shift/tab to edit Prefix and Suffix)
 - F1 for help website
 - Mouse-over "Help" button in upper-right for table of all shortcuts
- Interactive Help
 - "Did You Know...?" feature tutorials at launch
 - New "What's This?" help mode provides detailed contextual help on key features, including metadata on saved thumbnails
- Themes
 - Easy dark-vs-light mode, plus configurable skins / themes
 - Updated splashscreen provides version and load status

Laser Safety

- On-screen Status Indicators (virtual LEDs) reflect laser interlock and laser warning delay state on compatible devices
- Added second Laser Enable button in GUI masthead for visibility and access across all screens
- Configurable Laser Watchdog Delay on compatible devices
- Optional audio feedback on laser operation
- Ensure laser disabled when exiting Batch Collection
- Selectable laser power unit (% or mW) on ILC models and Expert mode

File Storage

- User is warned when renaming a thumbnail would overwrite existing files
- User is warned before erasing Clipboard
- Spectra can be saved in JCAMP-DX format
- Horizontal ROI
 - Processed data outside configured horizontal ROI saved as "NA" (raw spectra includes full detector area)
- Exports
 - Exports allow "non-collated" option grouped by spectrum type (processed, raw, dark etc)
 - Exports allow user to navigate to any folder when saving data
 - Configured Prefix and Suffix fields automatically applied to export filenames
 - Option to only export displayed traces
- Templating
 - Many additional {macro} types supported including time/date (see Manual)
 - On-screen filename Prefix and Suffix fields, as well as Notes metadata field, support {macro} template expansion

- User is warned when a "static" template lacking any {macro} fields for per-measurement customization would cause every measurement to be written to the same filename (useful in some streaming / watch-folder applications)
- Configurable thumbnail limit (default still 500)

Post-Processing

- ASTM-1840 Raman Wavenumber Shift Correction improvements
- Raman Intensity Correction (SRM) moved to Scope button-bar
- Interpolation moved to Scope button-bar

Product Support

- Support for 785XS
- Support for 1064XL
- Improved support for downloading XL "virtual EEPROMs" from cloud

Plugin Support

- Plugins are now accessed from the main scrolling right-hand Tools Palette (previously had been on Settings view)
- Plugin architecture has been streamlined to a more-Pythonic "functional" API using fewer classes and methods. All sample plugins have been updated to the new interface as demonstrations.
- Plugins now have access to all ENLIGHTEN internal state and functionality via a handle to the Controller
- A single Plugin can be auto-loaded at launch via command-line option
- New ComboBox field type
- Plugin fields can be marked for Expert mode
- Plugins can register themselves into Presets
- Plugin selection persisted across sessions
- Tensorflow support
- Sample Plugins
 - Raman.RamanPeaks shows expected peak locations for common compounds
 - Prod.EmissionLines shows expected atomic lines for noble gases

Manufacturing / Hardware Features

- Hardware View
 - EEPROM fields can be filtered via search bar

Cross-Platform

- Improved installation process (.dmg) for MacOS

- New "run from source" installation script for Linux (Ubuntu, Raspberry Pi)

Internal (Developer)

- Business Object constructors have been simplified by giving all objects a handle to the Controller
- Enforce DEBUG logging at launch, during spectrometer connection and 10sec thereafter

Appendix: ENLIGHTEN™ Changelog (3.2.36 - 4.1.0)

The following is extracted from the full ENLIGHTEN™ CHANGELOG.md available on GitHub at:

<https://github.com/WasatchPhotonics/ENLIGHTEN/blob/main/CHANGELOG.md>


4.0 GUI and Plugin Refactoring

- 2024-08-22 4.1.0
 - public release
- 2024-08-22 4.0.68
 - change StatusBar frame count from ProcessedReading to Reading.session_count
 - clarified BatchCollection "explain this" re: collection timeout
 - [#462] fix system clipboard copy after connecting to plugin
 - [#460] fix issues loading saved spectra
 - fix logging bug in MeasurementFactory
 - fixed missed calls to add_renamable from #452
 - fix bug graphing loaded .json measurements
 - enlighten_ubuntu_launcher.sh
 - fixed curl
 - added python3.7-dev
- 2024-07-19 4.0.67
 - [#452] don't prompt when re-labeling thumbnails
 - [#454] fix issue renaming .spc files
 - [#455] fix basename issue when multiple thumbnails have the same label
 - [#456] laser power unit combobox should disappear when leaving expert mode
 - [#450] indicate when laser in LaserWarningDelay state
- 2024-06-28 4.0.66
 - [#389] add Preset to metadata
 - [#406] handle whitespace in developer path
 - [#408] support --password cmd-line arg and ENLIGHTEN_PASSWORD environment variable
 - [#420] allow user to navigate to target directory when saving exports
 - [#426] support --max-thumbnails cmd-line arg and ENLIGHTEN_MAX_THUMBNAILS environment variable
 - [#428] verify before overwriting files
 - [#429] confirm when user configures "static" filename template
 - [#431] validate interpolation inputs
 - [#432] fixed scan averaging in BatchCollection
 - [#446] update theming docs
 - added hidden scipy imports for superman
- 2026-06-01 4.0.65
 - allow plugins to override Label metadata

- 2024-06-19 4.0.64
 - move pywin32 and spc_spectra to requirements.txt
 - try to remove console on Win11
 - add --noconsole to pyinstaller
 - add FakeOutputHandle to keep Tensorflow happy in absence of sys.stdout
 - struggle with Microsoft Defender Antivirus
 - played with Python 3.12, rolled-back to 3.11
 - accept bizarre 1min+ delay on first Tensorflow initialization post-install/compilation
- 2024-05-30 4.0.63
 - add default FWHM
 - plugins
 - Exporters/SaveAsAngstrom: modernized, simplified
 - Demo/EventResponse: modernized, simplified, renamed Demo/Buttons
 - Demo/LoadHook: modernized, simplified
 - Demo/MakeOddNegative: modernized, simplified
- 2024-05-22 4.0.62
 - XL
 - support horizontal ROI
 - logging
 - enforce DEBUG logging at launch
 - enforce DEBUG logging at spectrometer connection
 - relax to configured logging (typically INFO) 10sec after spectrometer connection
 - do NOT persist verbose logging through enlighten.ini
- 2024-05-14 4.0.61
 - changed laser status messages
- 2024-05-10 4.0.60
 - UI
 - display empirical assessment of laser firing status in Marquee
 - merged in mpylon's awesome EEPROMEditor search bar
 - Measurements
 - prompt user before overwriting existing files on export
 - add vertical ROI metadata
 - BatchCollection
 - always disable laser in BatchCollection.stop
 - added TakeOneFeature copy-ctor (fixed LaserMode.Spectrum)
 - plugins
 - fixed bug in base class
 - moved EmissionLines and RamanLines to combobox
 - XS
 - added microcontroller serial number
- 2024-04-05 4.0.59
 - fixed bug when exporting only displayed traces

- more-robust handling of corrupt EEPROMs
- master enlighten_ubuntu_launcher.sh within scripts
- plugins
 - moved PluginController to ctl
 - persist last-selected plugin
 - removed enlighten.Plugins.EnlightenApplicationInfoReal, EnlightenPlugin.EnlightenApplicationInfo
- 2024-03-25 4.0.58
 - restore saturation warning on XS
 - enable Marquee.error with "hazard" CSS
- 2024-03-25 4.0.57
 - fix check_refs in BoxcarFeature
 - support Expert mode in plugin fields
 - add Session Count to Measurement metadata
 - consistently handle plugin metadata
 - added PresetFeature.clear_all
- 2024-03-23 4.0.56
 - updated ColumnFileParser with wasatch.CSVLoader updates
 - changing Preset auto-clears Dark and Reference
 - fixed GUI.msgbox_with_checkbox
 - fixed SpectrometerApplicationState.check_refs
 - moved plugin persistence from "advanced_options" to "plugins"
 - SRM works properly in Presets
- 2024-03-22 4.0.55
 - added msgbox_with_lineedit
 - hopefully didn't break dialog response processing (narrator: it did)
 - R&D
 - skip ambient on ARM 1.0.2.9
 - ignore None responses in WasatchDevice funcs
- 2024-03-22 4.0.54
 - fixed GainDB persistence
 - fixed Reference visibility when Emission/Raman selected
 - flake8 fixes
 - moved GridFeature to ctl
- 2024-03-21 4.0.53
 - moved ReferenceFeature to ctl
 - fixed Dark/Reference on Settings
 - fixed Detector Temperature label bug
 - simplified designer.sh
 - added 'make flake8'
- 2024-03-19 4.0.52
 - started KIA conversion to ctl
 - fixed 2/3-byte FPGA_COMPILATION_OPTIONS
- 2024-03-19 4.0.51

- merged old Temperature Control panel into Detector Control widget (as Detector Temperature)
- added new Laser Temperature panel within Laser Control widget
- added Laser Temperature to StatusBar
- added EDC for IMX385
- LaserTemperatureFeature to ctl
 - added XS Laser TEC Mode
- added AmbientTemperatureFeature
- refactored Factory view a bit
- prompt before erasing Clipboard
- ctrl-5 for Factory view (authenticates if needed)
- only expose XS gain in Expert mode
- 2024-03-13 4.0.50
 - fixed changelog (2024)
 - SRM / Wavenumber button tweaks (also location)
 - ARM fixes regarding FPGAOptions and div-by-zero InGaAs correction
- 2024-03-11 4.0.49
 - DYK (Stats, StatsBuffer, EmissionLines, RamanLines)
 - added RamanIntensityCorrection button
 - added support for orange "tri-state" buttons
 - removed duplicate theme stylesheets
 - ROI editing tweaks
- 2024-03-07 4.0.48
 - temporarily interpolate clipboard measurements during export
 - display "sensor stabilizing" message on SiG
 - display wasatch.StatusMessage via Marquee, not Toast
- 2024-03-07 4.0.47
 - fixed export bug
- 2024-03-01 4.0.46
 - bumping version number for clean installer
- 2024-03-01 4.0.45
 - Multispec fix
 - Tensorflow Win11 fix
 - export ROI fix
 - flipped (corrected) meaning of "collated_export" (still defaults to "original behavior")
 - restore installer COMPRESSION option
- 2024-02-27 4.0.44
 - deprecated Miniconda (python venv FTW)
 - built installer on Win11
 - ASTM tweak
 - DYK tweaks
 - convert timestamps to strings in JSON
 - fix pixel axis with ROI

- disable KIA tips
- 2024-02-20 4.0.43
 - fixed XL temperature
 - fixed Peak Sharpening
- 2024-02-19 4.0.42
 - GUI tweaks
 - added DidYouKnowFeature
 - added session close button for KIA
 - fleshed-out WhatsThis docs
 - theme fixes
 - moved sounds checkbox to Settings view
 - added laser_steady.wav
 - tightened checks around StatusIndicators
- 2024-02-14 4.0.41 
 - RamanShiftCorrection: better support for ROI, tweaked peak windows
 - updated to Tensorflow 2.13.1
- 2024-02-09 4.0.40
 - Configuration cleanup
 - persist scans to average in enlighten.ini by serial number
 - move PluginController.autoload to PluginController QTimer
- 2024-02-07 4.0.39
 - fixed ROI bug in Measurement.save_excel_file
 - properly apply ROI for loaded export and columnar csv files
- 2024-02-06 4.0.38
 - fixed load from export
 - fixed graph traces on loaded exports
 - fixed common.msgbox bug
 - moved to ctl
 - ColumnFileParser
 - DashFileParser
 - ExportFileParser
 - SPCFileParser
 - TextFileParser
- 2024-02-05 4.0.37
 - pylint
 - docs
 - fixed LaserControlFeature restriction bug
 - fixed Measurements export ROI bug
- 2024-01-30 4.0.36
 - Raman
 - ensured all Raman features available in Expert mode
 - renamed RamanModeFeature to AutoRaman
 - moving Raman Intensity Correction to its own widget visible in Raman view

- extracted LaserWatchdogFeature
 - don't force Laser Watchdog unless XS SML w/o TEC (but support for any XS)
- GUI
 - Guide tweaks
 - common.msgbox tweaks
- cleanup
 - moved Filters and Authentication to ui
 - moved SpectrometerApplicationState, AccessoryControlFeature, DetectorTemperatureFeature, ExternalTriggerFeature, GainDBFeature, HighGainModeFeature, IntegrationTimeFeature, LaserControlFeature, LaserTemperatureFeature, MultiPos, RegionControlFeature and Spectrometer to device
 - renamed spectra_processes to post_processing
 - moved DarkFeature, ReferenceFeature, RamanMode to post_processing
 - renamed pluginExamples to plugins
 - renamed sfu to cfu
- moved to ctl
 - HighGainModeFeature
 - ManufacturingFeature
 - CloudManager
 - RichardsonLucy
 - Sounds
- 2024-01-26 4.0.35
 - RamanModeFeature
 - prompt for confirmation
 - release restriction on LaserControlFeature when invisible
 - better handling of case where laser already firing when Raman Mode engaged
 - BatchCollection
 - bugfix (typo)
 - disable Guide during BatchCollection
 - enforce 2sec gap between subsequent Batches to allow manual "Stop" and stabilization
 - reset TakeOneFeature before calling TakeOne.complete observers
 - added --start-batch cmd-line
 - Plugins
 - simplified --plugin cmd-line arg
 - made RamanLines auto-enable
 - HelpFeature: added shortcuts to tooltip
 - LoggingFeature: utf-8 bugfix
- 2024-01-24 4.0.34
 - Refactored RamanModeFeature, VCRControls and TakeOneFeature around TakeOneRequest

- moved to ctl
 - BatchCollection
 - RamanModeFeature
 - TakeOneFeature
 - VCRControls
- 2024-01-15 4.0.33
 - fix load bug
 - simplify common.msgbox api
 - simplify thumbnail removal prompt
 - correctly associate thumbnails with generating plugin
 - give Graphs names for tracking
 - make utf8 settings more explicit
 - default logfile to truncate
- 2024-01-12 4.0.32
 - GUI
 - StatusBar only shows battery when appropriate
 - added common.msgbox(..., detail)
 - added GUI.msgbox_with_textedit
 - plugins
 - renamed logfile to plugin.log
 - truncate plugin logfile on load
 - updated peaks in RamanLines and ASTM-1840 configuration
- 2024-01-11 4.0.31
 - axis corner-cases
 - replaced InterpolatedProcessedReading with ProcessedReading.cropped and .interpolated
 - prevent export from trying to mix different x-axes from the same serial number
 - GUI
 - allow visualizing SRM calibration
 - styled common.msgbox
 - expand ScrollStealFilter to include focused widgets
 - fixed bug graphing pixel axis
 - moved to ctl
 - Authentication
 - EEPROMWriter
 - MeasurementFactory
 - scripts
 - retain more metadata in split-spectra
- 2023-12-14 4.0.30
 - GUI
 - added battery to StatusBar
 - generalized PluginController.display_plugin_warning into GUI.msgbox_with_checkbox

- simplified font sizes
 - added X series png's
 - fixed Hardware tab device image
 - optionally Export only displayed traces
 - moved to ctl
 - AbsorbanceFeature
 - Clipboard
 - Configuration
 - DetectorTemperatureFeature
 - EEPROMEditor
 - ExternalTriggerFeature
 - FileManager
 - FocusListener
 - GuideFeature
 - HardwareFileOutputManager
 - LoggingFeature
 - Marquee
 - ModelInfo
 - Multispec
 - ResourceMonitorFeature
 - SaveOptions
 - StatusBarFeature
 - ThumbnailWidget
 - TransmissionFeature
 - scripts
 - refreshed scripts/split-spectra
- 2023-12-13 4.0.29
 - stopped clearing dark/refs when changing interpolation
 - RamanShiftCorrection tweaks (fixed width unit, obey ROI, show nm hint)
 - refactoring
 - moved DarkFeature to ctl
- 2023-12-11 4.0.28
 - platforms
 - updating Raspberry Pi build docs
 - support both PySide6 and PySide2
 - GUI
 - laser button tooltip fix
 - fixed markers
 - added convenience laser button where permanently visible
 - added interpolation button above graph
 - restored original ordering of (lock, fullwidth) buttons
 - fixed multiple Factory items on multiple logins
- 2023-12-08 4.0.27
 - locked to Python 3.11 for Tensorflow

- default to maximized
- restore previous integration time up to 5sec (was 1sec)
- plugins
 - all field types support callbacks
- docs
 - updated Windows docs for 10/11
 - updated MacOS docs for ARM
- 2023-11-10 4.0.26
 - auto-resize pandas table columns
 - fix ctl issues in Measurement*
 - move jcamp from pip to PYTHONPATH until feature branch released
- 2023-11-06 4.0.25
 - added PresetFeature
 - plugins
 - hide plugin field frame if no fields
 - deprecate dependencies
 - support events in functional plugins
 - increase default plugin Marquee period
 - provide self.marquee_message for functional plugins
 - add combobox datatype for plugin fields
 - allow plugin fields to have stylesheets
 - auto-resize pandas table columns
 - file_io
 - persist custom spectra save location
 - save as JCAMP-DX
 - support loading columnar CSV files which don't happen to include 'Integration Time'
 - refactoring
 - moved RamanIntensityCorrection, ScanAveragingFeature, BaselineCorrection, BoxcarFeature, IntegrationTimeFeature to ctl
 - docs and cleanup around ScrollStealFilter and MouseWheelFilter
 - simplified Configuration.save
- 2023-10-09 4.0.24
 - Prod.BurnIn temperature fix
 - added run_from_root
 - restored hide_from_console (hide-early doesn't work on W11, hide-late doesn't on W10)
- 2023-10-06 4.0.23
 - temporarily changed PluginController autoloader from "ctor + 1sec" to "first spectrometer connection" (later rolled-back, as broke ID plugin)
- 2023-10-06 4.0.22
 - bundle libusb0.dll to Windows install folder :-)
- 2023-10-06 4.0.21
 - automatically retain previous 2MB of existing enlighten.log to aid in fault analysis

- 2023-09-30 4.0.20
 - updated dependency versions to resolve tensorflow import issue
 - Python 3.10 -> 3.11
 - PySide 2 -> 6
 - pyinstaller 4 -> 5
 - numpy -> 1.26 (1.24.3?)
 - tensorflow <- 2.13
 - let bootstrap pass compression to InnoSetup (default lzma/max)
 - change pyinstaller --noconsole to --hide-console to fix pexpect issue
- 2023-09-29 4.0.19
 - applog tweak
- 2023-09-29 4.0.18
 - fixed plugin outputs/table bug
 - added ScanAveragingFeature.set_locked
 - changed --headless to --window-state
 - BurnIn tweak
- 2023-09-20 4.0.17
 - support --log-append
 - moved run-burning.bat to plugins/Prod
 - improved Prod.BurnIn logging
- 2023-09-20 4.0.16
 - bugfixes
 - RamanIntensityCorrection works in Expert Mode
 - use EEPROM FWHM for RamanShiftCorrection
 - UI
 - splash tweak
 - msgbox for improved error reporting (e.g. 0 max integration time)
 - portability
 - fixed rebuild_resources.sh for MacOS
 - updated Mac .dmg installer process
 - production
 - added BurnIn plugin
- 2023-09-01 4.0.15
 - enforce eeprom.max_integration_time_ms on XS
 - fixed DetectorTemperatureFeature on ambient detectors
 - plugins
 - added --plugin command-line option
 - added Demo.LockSettings
 - XS-Series
 - enforce eeprom.max_integration_time_ms
 - StatusIndicators show laser unavailable when battery low
 - refactoring
 - moved BatteryFeature, StatusIndicators to ctl
- 2023-07-19 4.0.14

- extend boxcar to left and right edges of spectra
- clear scan average count upon reset
- 2023-07-07 4.0.13
 - spectrometer control
 - support selectable mW/PWM for ILC and Expert Mode
 - data management
 - added "use filename as label" checkbox
 - template improvements
 - added filename template
 - macros now work in export filename components
 - support Reading attributes
 - user-configurable floating-point precision
 - added date/time components ({YYYY}, {hh} etc)
 - added file_timestamp, integration_time_sec macros
 - persistence
 - enlighten.ini allows %
 - persist gain_db by serial number
 - GUI updates
 - added WhatsThis button next to Help
 - changed Measurement ToolTip to WhatsThis
 - added min/avg/max to thumbnail WhatsThis
 - thumbnail label improvements (smaller font, easier editing)
 - keyboard shortcuts
 - added ctrl-G/N/T per <https://www.goodreads.com/quotes/9934315-it-is-a-curious-fact-and-one-to-which-no>
 - added F1 for help
 - splash screen messages
 - plugin improvements
 - updated Savitzky-Golay
 - updated StatsBuffer
 - updated Stats
 - updated docs
 - fixed plugin_name in metadata
 - refactoring
 - renamed VignetteROIFeature to HorizROIFeature
 - updated LaserControlFeature to ctl
 - technique-vs-view fixes
 - vignette-vs-crop fixes
 - static analysis support
 - devtools
 - pyusb-virtSpec support
- 2023-06-14 4.0.12
 - persistence

- fixed Configuration defaults
 - updated Write EEPROM tooltip
 - file_io
 - warn when saving but no formats selected
 - fixed load-then-export use-case
 - always save all rows of column-ordered CSV spectra regardless of ROI
 - output "NA" rather than 0 for values omitted due to ROI or spectral range
 - use prefix and suffix in export filenames
 - stop forcing loaded data to lowercase
 - improved data consistency in loaded CSVs
 - plugins
 - moved from Settings View to Control Palette
 - moved RamanLines to Raman
 - simplified Prod.EmissionLines
- 2023-05-15 4.0.11
 - added collated export
 - invert x-axis for 1064XL
- 2023-05-10 4.0.10
 - "fixed" Marquee
 - added scripts/transpose.py
 - fixed save on non-Raman
 - allow LocalBaseline plugin to save metadata
 - persist cloud permission
 - fixed Raman Intensity Correction in 1064XL
- 2023-05-03 4.0.9
 - fixed Save and Export
- 2023-05-02 4.0.8
 - added LocalBaseline plugin
 - refactored Business Objects package structure
 - pyinstaller tweak (excluded _bootlocale)
- 2023-03-27 4.0.7
 - persist PluginController marquee while satisfying prerequisites
 - clamped Python to 3.10 to resolve build issues
 - removed matplotlib dependency from non-Mac platforms
 - added cyclohexane peaks
 - fixed PYTHONPATH issue affecting installer
- 2023-03-17 4.0.6
 - Added Plugin Quick Start documentation
 - Added Worek AChE Timeseries Plugin
 - Ctrl-E to rename
 - Fix bugs causing CSV files to fail to load
 - Added Hello Graph Demo Plugin
- 2023-02-08 4.0.5
 - watchdog toggle via checkbox

- minimum watchdog of 3 sec
- 2023-01-31 4.0.4
 - fixed logo in light mode
 - fixed installer bug?
- 2023-01-13 4.0.3
 - added dark-mode button
 - hopefully improved laser watchdog handling on older (unsupported) models
- 2023-01-03 4.0.2
 - persist initial spectrometer color assignment across disconnects within a session
 - added thumbnail tooltip
- 2022-12-19 4.0.1
 - improved plugin exception-handling
 - improved plugin dependency messaging
 - improved Marquee 'immediate' display
- 2022-12-09 4.0.0
 - better widget scaling at lower resolutions
 - refactored Techniques, Modes and Views
 - refactored GridFeature and VignetteROIFeature
 - refactored LoggingFeature
 - reflect laserWatchdog state in GUI
 - fixed bug when multiple thumbnails have same label

3.2 Open-Source

- 2022-12-08 3.2.42
 - expose save options to plugins
 - plugins can programmatically change input fields
 - plugins can save in same date dir as regular save
 - fix graphing issues for plugins
 - add measurement metadata so plugin thumbnail overlays show on other graph
- 2022-12-06 3.2.41
 - adding widget for laser watchdog
- 2022-11-29 3.2.40
 - dark-correct reference spectrum after-the-fact if collected in that order
 - don't apply boxcar to dark or reference when in reference technique
- 2022-11-28 3.2.39
 - i18n
 - add fix for international computers for plugins creation on install
- 2022-11-15 3.2.38
 - add battery data to hardware strip charts
 - add string formatting to note, prefix, and suffix
- 2022-10-17 3.2.37
 - updates to regression tests

- 2022-10-12 3.2.36
 - merged pluginWarning / ROI fixes

Appendix: Wasatch.PY Changelog (2.1.21 - 2.2.0)

The following is extracted from the full Wasatch.PY CHANGELOG.md available on GitHub at:

<https://github.com/WasatchPhotonics/Wasatch.PY/blob/master/CHANGELOG.md>

- 2024-08-22 2.2.0
 - ENLIGHTEN 4.1.0 release
- 2024-06-27 2.1.73
 - update demo.py to support multiple spectrometers
 - add optional callback parameter to WrapperWorker for Inversion of Control (IoC)
- 2024-06-27 2.1.72
 - temporarily disable GET_LASER_ENABLED on XS
- 2024-06-19 2.1.71
 - disable applog to stdout if stdout missing
- 2024-06-12 2.1.70
 - force laser to disable on perc/PWM zero
- 2024-05-30 2.1.69
 - added avg_resolution to editable EEPROM fields
- 2024-05-24 2.1.68
 - added horizontal ROI to XL virtual EEPROM
 - Raspberry Pi testing (no change)
- 2024-05-22 2.1.67
 - relax thread logging 10sec after successful connection
 - display connected spectrometer properties in demo.py
- 2024-05-07 2.1.66
 - gave TakeOneRequest a copy-ctor
- 2024-05-02 2.1.65
 - updated README links to work on pypi servers
- 2024-05-02 2.1.64
 - added wasatch.FirmwareRequirements
 - added wasatch.PollStatus
 - added FID.get_poll_status() and .is_sensor_stable()
 - added util.vercmp()
 - added microcontroller serial number
- 2024-03-23 2.1.63
 - CSVLoader
 - store x-axes in ProcessedReading rather than metadata
 - support "intensity" header
 - simplified
- 2024-03-22 2.1.62
 - ignore 3 consecutive USB errors when reading spectra (R&D)
 - ignore None responses in WasatchDevice func/attr calls
 - skip ambient temperature readout on ARM 1.0.2.9
- 2024-03-19 2.1.61

- o handle either 3-byte or 2-byte FPGA_COMPILATION_OPTIONS response
- 2024-03-18 2.1.60
 - o added get_ambient_temperature_degC_arm
- 2024-03-15 2.1.59
 - o changed XS get_detector_gain to big-endian
 - o added clarity around XS laser_tec_mode -vs- enabled
- 2024-03-15 2.1.58
 - o added set/get_laser_tec_mode
 - o sanity-check XS laser TEC setpoint to 700-900 range
- 2024-03-14 2.1.57
 - o added laser temperature calibration for XS
 - o added EDC for XS
 - o clarified FPGAOptions as 24-bit (but ignore for XS)
- 2024-03-13 2.1.56
 - o double-check to prevent div-by-zero in InGaAs gain correction
 - o decline to parse FPGAOptions for ARM (leave at initialized values)
- 2024-03-13 2.1.55
 - o don't use FPGAOptions on ARM
 - o return "Sensor Stabilizing" message on ARM when intTime/gain changes
- 2024-03-07 2.1.54
 - o return "sensor stabilizing" message
- 2024-02-28 2.1.53
 - o enriched ProcessedReading.dump
- 2024-02-28 2.1.52
 - o fixed ProcessedReading bug
- 2024-02-27 2.1.51
 - o added ProcessedReading.get_pixel_axis
- 2024-02-20 2.1.50
 - o fixed AndorDevice temperature
- 2024-02-19 2.1.49
 - o return state.laser_enabled from FID.get_laser_enabled on SiG until FW fixed
- 2024-02-14 2.1.48
 - o ProcessedReading: force to np.array on load_from_dict, clear .interpolated on set_processed
- 2024-02-07 2.1.47
 - o ProcessedReading always takes wavelengths/wavenumbers from settings if empty and available
- 2024-02-06 2.1.46
 - o removed enlighten.SaveOptions(?) from wasatch.CSVLoader(?!?)
 - o made ProcessedReading smarter about loading wavelengths/wavenumbers from SpectrometerSettings
- 2024-02-05 2.1.45
 - o clarified scan averaging a bit
 - o added TakeOneRequest.scans_to_average
- 2024-01-30 2.1.44

- o Auto-Raman fix: store laser-related attributes in Reading before auto-disabling laser
 - o notation updates
- 2024-01-26 2.1.43
 - o utils bugfix
- 2024-01-25 2.1.42
 - o TakeOneRequest.take_dark + .enable_laser_before now implies laser disable + 1sec delay
- 2024-01-23 2.1.41
 - o moved to TakeOneRequest
 - o deprecated acquisition_laser_trigger_enable, self.acquisition_take_dark_enable and acquisition_laser_trigger_delay_ms
- 2024-01-15 2.1.40
 - o better UTF-8 support
 - o support cropped/interpolated ProcessedReading in dict conversions
- 2024-01-11 2.1.39
 - o refactored ProcessedReading with new .cropped and .interpolated nodes
 - o fixed small PWM bug
- 2023-12-18 2.1.38
 - o fixed iterable bug in ProcessedReading dict ctor
- 2023-12-14 2.1.37
 - o ProcessedReading tweak
- 2023-12-07 2.1.36
 - o add throwaways around IMX sensor commands
 - o treat flat IMX spectra as keepalives
- 2023-11-01 2.1.35
 - o auto-generate raman_intensity_calibration_order in AndorDevice virtual EEPROM if missing
 - o added ProcessedReading.to_dict()
- 2023-09-20 2.1.34
 - o added applog(append)
- 2023-09-20 2.1.33
 - o enforce FPGA version number length
 - o log USB results in hex
 - o check for "bad data" in EEPROM
 - o XS-Series
 - added set_laser_power_attenuator
 - added support for startup laser TEC setpoint
 - read battery @ 0.2Hz
- 2023-07-10 2.1.32
 - o offset getters now return SpectrometerResponse
 - o fixed display of getters returning SpectrometerResponse in WasatchShell
- 2023-07-04 2.1.31
 - o clarify methods returning bool vs SpectrometerResponse
 - o remove return type annotations
 - o s/vignette/crop/

- 2023-06-13 2.1.30
 - invert x-axis for 1064XL
 - don't enable TEC with default / invalid TEC calibration
 - restore ability for CSVLoader.parse_metadata to support single-measurement CSVs where metadata fields are expected to be scalars rather than lists
- 2023-05-10 2.1.29
 - support for 1064XL
- 2023-05-03 2.1.28
 - better handling of rare null-buffer EEPROM readout
- 2023-01-13 2.1.27
 - is_andor fix
 - laser watchdog fix (disable on older EEPROMs)
- 2022-12-14 2.1.26
 - default InGaAs detectors to high-gain mode
 - forget disconnected USB devices
- 2022-12-13 2.1.25
 - disable laserEnable verification loop for Series-XS
- 2022-12-07 2.1.24
 - update EEPROM to Rev15
 - default laser power calibration to high-resolution
- 2022-10-26 2.1.23
 - fix connect bug
 - added scripts/dark-vs-integ.py
- 2022-10-17 2.1.22
 - updates to MockUSBDevice to support ENLIGHTEN regression tests
- 2022-10-12 2.1.21
 - add is_sml and is_mml functions
 - fix floating-point array index lookup